

Sanitize

Complete Track and Index Listing

<u>Filename</u>	<u>Description</u>	<u>Duration</u>	<u>Library</u>	<u>Category</u>	<u>Sample Rate</u>	<u>Bit Depth</u>
SAN001	Purell bottle open and close, handling		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN002	Sanitizer foot operated stand with bottle		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN003	Sanitizer foot operated stand without bottle		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN004	Pump action sanitizer bottle empty		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN005	Pump action sanitizer bottle full		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN006	Sanitize hand rubbing		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN007	Spray bottle small		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN008	Spray bottle large		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN009	Disposable material cloth handle and tear		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN010	Disposable paper towel handle and tear		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN011	Rubber gloves		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN012	Wet wipes removal and handling		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN013	Tissue box		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN014	Face mask handling, material		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN015	Face shield handling		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN016	Infrared thermometer		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN017	Handwashing ceramic basin with water running		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN018	Handwashing no water		TWRSFX021 - Sanitize	Foley	96Khz	24Bit
SAN019	Handwashing steel basin with water running		TWRSFX021 - Sanitize	Foley	96Khz	24Bit



THEWORKROOM
Audio Post Production